

For Research Use Only

HIC2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00206-2

Capture Antibody Information

Catalog Number:
83162-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230177D1
Reactivity:
Human
GenBank:
BC094787
Immunogen Catalog Number:
Ag18789

Conjugate:
Unconjugated
Full name:
hypermethylated in cancer 2
Gene ID:
23119

Detection Antibody Information

Catalog Number:
83162-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230177H1
Reactivity:
Human
GenBank:
BC094787
Immunogen Catalog Number:
Ag18789

Conjugate:
Unconjugated
Full name:
hypermethylated in cancer 2
Gene ID:
23119

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00206-2 targets HIC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HIC2 Recombinant antibody, PBS Only (Capture) 83162-3-PBS (230177D1). 100 µg. Concentration 1 mg/ml.

Detection antibody: HIC2 Recombinant antibody, PBS Only (Detector) 83162-4-PBS (230177H1). 100 µg. Concentration 1 mg/ml.

Alternative HIC2 matched antibody pairs: MP00206-1, MP00206-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00206-2, HIC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83162-3-PBS. Detection antibody: 83162-4-PBS. Standard: Ag18789. Range: 0.625-20 ng/mL