

For Research Use Only

HIC2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00206-3

Capture Antibody Information

Catalog Number: 83162-6-PBS	Clone ID: 230177D8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: hypermethylated in cancer 2
Isotype: IgG	GenBank: BC094787	Gene ID: 23119
Purification Method: Protein A purification	Immunogen Catalog Number: Ag18789	

Detection Antibody Information

Catalog Number: 83162-2-PBS	Clone ID: 230177G7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: hypermethylated in cancer 2
Isotype: IgG	GenBank: BC094787	Gene ID: 23119
Purification Method: Protein A purification	Immunogen Catalog Number: Ag18789	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00206-3 targets HIC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HIC2 Recombinant antibody, PBS Only (Capture) 83162-6-PBS (230177D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: HIC2 Recombinant antibody, PBS Only (Detector) 83162-2-PBS (230177G7). 100 µg. Concentration 1 mg/ml.

Alternative HIC2 matched antibody pairs: MP00206-1, MP00206-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00206-3, HIC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83162-6-PBS. Detection antibody: 83162-2-PBS. Standard: Ag18789. Range: 0.625-20 ng/mL