

For Research Use Only



TMEM65 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00209-1

Capture Antibody Information

Catalog Number: 83219-2-PBS	Clone ID: 230512E3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: transmembrane protein 65
Isotype: IgG	GenBank: BC032396	Gene ID: 157378
Purification Method: Protein A purification	Immunogen Catalog Number: Ag15509	

Detection Antibody Information

Catalog Number: 83219-1-PBS	Clone ID: 230512A3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: transmembrane protein 65
Isotype: IgG	GenBank: BC032396	Gene ID: 157378
Purification Method: Protein A purification	Immunogen Catalog Number: Ag15509	

Applications

Tested Applications: Cytometric bead array	Range: 6.25-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00209-1 targets TMEM65 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMEM65 Recombinant antibody, PBS Only (Capture) 83219-2-PBS (230512E3). 100 μ g. Concentration 1 mg/ml.

Detection antibody: TMEM65 Recombinant antibody, PBS Only (Detector) 83219-1-PBS (230512A3). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00209-1, TMEM65 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83219-2-PBS. Detection antibody: 83219-1-PBS. Standard: Ag15509. Range: 6.25-100 ng/mL.