

RILPL1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00211-1

Capture Antibody Information

Catalog Number:
83220-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230515D7
Reactivity:
Human
GenBank:
BC080626
Immunogen Catalog Number:
Ag10192

Conjugate:
Unconjugated
Full name:
Rab interacting lysosomal protein-like 1
Gene ID:
353116

Detection Antibody Information

Catalog Number:
83220-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230515G8
Reactivity:
Human
GenBank:
BC080626
Immunogen Catalog Number:
Ag10192

Conjugate:
Unconjugated
Full name:
Rab interacting lysosomal protein-like 1
Gene ID:
353116

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00211-1 targets RILPL1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RILPL1 Recombinant antibody, PBS Only (Capture) 83220-5-PBS (230515D7). 100 µg. Concentration 1 mg/mL.

Detection antibody: RILPL1 Recombinant antibody, PBS Only (Detector) 83220-2-PBS (230515G8). 100 µg. Concentration 1 mg/mL.

Alternative RILPL1 matched antibody pairs: MP00211-2, MP00211-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00211-1, RILPL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83220-5-PBS. Detection antibody: 83220-2-PBS. Standard: Ag10192. Range: 0.625-20 ng/mL