

For Research Use Only

MAP7D1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00212-3

Capture Antibody Information

Catalog Number:
83173-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230424G2
Reactivity:
Human
GenBank:
BC003083
Immunogen Catalog Number:
Ag14385

Conjugate:
Unconjugated
Full name:
MAP7 domain containing 1
Gene ID:
55700

Detection Antibody Information

Catalog Number:
83173-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230424D8
Reactivity:
Human
GenBank:
BC003083
Immunogen Catalog Number:
Ag14385

Conjugate:
Unconjugated
Full name:
MAP7 domain containing 1
Gene ID:
55700

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00212-3 targets MAP7D1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MAP7D1 Recombinant antibody, PBS Only (Capture) 83173-4-PBS (230424G2). 100 µg. Concentration 1 mg/ml.

Detection antibody: MAP7D1 Recombinant antibody, PBS Only (Capture/Detector) 83173-3-PBS (230424D8). 100 µg. Concentration 1 mg/ml.

Alternative MAP7D1 matched antibody pairs: MP00212-1, MP00212-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00212-3, MAP7D1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83173-4-PBS. Detection antibody: 83173-3-PBS. Standard: Ag14385. Range: 0.625-20 ng/mL