

For Research Use Only

Rat KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00213-2

Capture Antibody Information

Catalog Number:
83221-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230528C5

Reactivity:
rat, mouse, human

GenBank:
NM_173149.2

Conjugate:
Unconjugated

Full name:
hepatitis A virus cellular receptor 1

Gene ID:
286934

Detection Antibody Information

Catalog Number:
83221-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230528G9

Reactivity:
rat

GenBank:
NM_173149.2

Conjugate:
Unconjugated

Full name:
hepatitis A virus cellular receptor 1

Gene ID:
286934

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00213-2 targets KIM-1/HAVCR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat KIM-1/HAVCR1 Recombinant antibody, PBS Only (Capture) 83221-2-PBS (230528C5). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Rat KIM-1/HAVCR1 Recombinant antibody, PBS Only (Detector) 83221-4-PBS (230528G9). 100 μ g. Concentration 1 mg/ml.

Alternative KIM-1/HAVCR1 matched antibody pairs: MP00116-1, MP00116-2, MP00116-3, MP00213-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00213-2, Rat KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83221-2-PBS. Detection antibody: 83221-4-PBS. Standard: Eg0597. Range: 0.625-20 ng/mL.