

For Research Use Only

MKKNK2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00216-1

Capture Antibody Information

Catalog Number:
83235-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240045D11

Reactivity:
human

GenBank:
BC018345

Immunogen Catalog Number:
Ag10546

Conjugate:
Unconjugated

Full name:
MAP kinase interacting
serine/threonine kinase 2

Gene ID:
2872

Detection Antibody Information

Catalog Number:
83235-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240045B9

Reactivity:
human

GenBank:
BC018345

Immunogen Catalog Number:
Ag10546

Conjugate:
Unconjugated

Full name:
MAP kinase interacting
serine/threonine kinase 2

Gene ID:
2872

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00216-1 targets MKKNK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MKKNK2 Recombinant antibody, PBS Only (Capture/Detector) 83235-4-PBS (240045D11). 100 μ g. Concentration 1 mg/ml.

Detection antibody: MKKNK2 Recombinant antibody, PBS Only (Capture/Detector) 83235-3-PBS (240045B9). 100 μ g. Concentration 1 mg/ml.

Alternative MKKNK2 matched antibody pairs: MP00216-2, MP00216-3, MP00216-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00216-1, MKNK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83235-4-PBS. Detection antibody: 83235-3-PBS. Standard: Ag10546. Range: 1.25-40 ng/mL