

For Research Use Only

MKKNK2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00216-3

Capture Antibody Information

Catalog Number:
83235-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240045B9
Reactivity:
human
GenBank:
BC018345
Immunogen Catalog Number:
Ag10546

Conjugate:
Unconjugated
Full name:
MAP kinase interacting serine/threonine kinase 2
Gene ID:
2872

Detection Antibody Information

Catalog Number:
83235-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240045C7
Reactivity:
human
GenBank:
BC018345
Immunogen Catalog Number:
Ag10546

Conjugate:
Unconjugated
Full name:
MAP kinase interacting serine/threonine kinase 2
Gene ID:
2872

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00216-3 targets MKKNK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MKKNK2 Recombinant antibody, PBS Only (Capture/Detector) 83235-3-PBS (240045B9). 100 μ g. Concentration 1 mg/ml.

Detection antibody: MKKNK2 Recombinant antibody, PBS Only (Detector) 83235-2-PBS (240045C7). 100 μ g. Concentration 1 mg/ml.

Alternative MKKNK2 matched antibody pairs: MP00216-1, MP00216-2, MP00216-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00216-3, MKNK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83235-3-PBS. Detection antibody: 83235-2-PBS. Standard: Ag10546. Range: 1.25-40 ng/mL