

For Research Use Only

# ABHD6 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00225-1

## Capture Antibody Information

**Catalog Number:**  
83249-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230472C1  
**Reactivity:**  
Human  
**GenBank:**  
BC001698  
**Immunogen Catalog Number:**  
Ag13968

**Conjugate:**  
Unconjugated  
**Full name:**  
abhydrolase domain containing 6  
**Gene ID:**  
57406

## Detection Antibody Information

**Catalog Number:**  
83249-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230472F9  
**Reactivity:**  
Human  
**GenBank:**  
BC001698  
**Immunogen Catalog Number:**  
Ag13968

**Conjugate:**  
Unconjugated  
**Full name:**  
abhydrolase domain containing 6  
**Gene ID:**  
57406

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
3.125-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00225-1 targets ABHD6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ABHD6 Recombinant antibody, PBS Only (Capture) 83249-1-PBS (230472C1). 100 µg. Concentration 1 mg/ml.

Detection antibody: ABHD6 Recombinant antibody, PBS Only (Detector) 83249-2-PBS (230472F9). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00225-1, ABHD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83249-1-PBS. Detection antibody: 83249-2-PBS. Standard: Ag13968. Range: 3.125-100 ng/mL.