

TRIM54 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00236-1

Capture Antibody Information

Catalog Number:
83175-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230391A4
Reactivity:
Human
GenBank:
BC141807
Immunogen Catalog Number:
Ag15229

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 54
Gene ID:
57159

Detection Antibody Information

Catalog Number:
83175-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230391F10
Reactivity:
Human
GenBank:
BC141807
Immunogen Catalog Number:
Ag15229

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 54
Gene ID:
57159

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00236-1 targets TRIM54 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM54 Recombinant antibody, PBS Only (Capture) 83175-1-PBS (230391A4). 100 µg. Concentration 1 mg/mL.

Detection antibody: TRIM54 Recombinant antibody, PBS Only (Detector) 83175-5-PBS (230391F10). 100 µg. Concentration 1 mg/mL.

Alternative TRIM54 matched antibody pairs: MP00236-2, MP00236-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00236-1, TRIM54 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83175-1-PBS. Detection antibody: 83175-5-PBS. Standard: Ag15229. Range: 0.625-20 ng/mL