

For Research Use Only

PLK1-phospho Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00238-1

Capture Antibody Information

Catalog Number:
83260-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240092D5

Reactivity:
Human

GenBank:
NM_005030

Conjugate:
Unconjugated

Full name:
polo-like kinase 1 (Drosophila)

Gene ID:
5347

Detection Antibody Information

Catalog Number:
83260-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240092A12

Reactivity:
Human

GenBank:
NM_005030

Conjugate:
Unconjugated

Full name:
polo-like kinase 1 (Drosophila)

Gene ID:
5347

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-10 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00238-1 targets PLK1-phospho in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PLK1-phospho Recombinant antibody, PBS Only (Capture) 83260-2-PBS (240092D5). 100 μ g. Concentration 1 mg/ml.

Detection antibody: PLK1-phospho Recombinant antibody, PBS Only (Detector) 83260-3-PBS (240092A12). 100 μ g. Concentration 1 mg/ml.

Alternative PLK1-phospho matched antibody pairs: MP00238-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00238-1, PLK1-phospho Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83260-2-PBS. Detection antibody: 83260-3-PBS. Standard: SY00401. Range: 0.156-10 ng/mL.