

For Research Use Only

MS4A12 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00244-2

Capture Antibody Information

Catalog Number:
83261-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240105D5
Reactivity:
Human
GenBank:
BC029793
Immunogen Catalog Number:
Ag3968

Conjugate:
Unconjugated
Full name:
membrane-spanning 4-domains, subfamily A, member 12
Gene ID:
54860

Detection Antibody Information

Catalog Number:
83261-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240105A11
Reactivity:
Human
GenBank:
BC029793
Immunogen Catalog Number:
Ag3968

Conjugate:
Unconjugated
Full name:
membrane-spanning 4-domains, subfamily A, member 12
Gene ID:
54860

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00244-2 targets MS4A12 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MS4A12 Recombinant antibody, PBS Only (Capture) 83261-2-PBS (240105D5). 100 µg. Concentration 1 mg/mL.

Detection antibody: MS4A12 Recombinant antibody, PBS Only (Detector) 83261-1-PBS (240105A11). 100 µg. Concentration 1 mg/mL.

Alternative MS4A12 matched antibody pairs: MP00244-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00244-2, MS4A12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83261-2-PBS. Detection antibody: 83261-1-PBS. Standard: Ag3968. Range: 0.625-20 ng/mL.