

For Research Use Only

TIMP1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00245-2

Capture Antibody Information

Catalog Number:
83151-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230329C5

Reactivity:
Human

GenBank:
BC000866

Conjugate:
Unconjugated

Full name:
TIMP metalloproteinase inhibitor 1

Gene ID:
7076

Affinity:
 $K_D = 6.82 \times 10^{-11}$
 $K_{off} = 1.35 \times 10^{-4}$
 $K_{on} = 1.98 \times 10^6$

Detection Antibody Information

Catalog Number:
83151-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230329F11

Reactivity:
Human

GenBank:
BC000866

Conjugate:
Unconjugated

Full name:
TIMP metalloproteinase inhibitor 1

Gene ID:
7076

Affinity:
 $K_D = 1.75 \times 10^{-11}$
 $K_{off} = 5.52 \times 10^{-6}$
 $K_{on} = 3.15 \times 10^5$

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00245-2 targets TIMP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TIMP1 Recombinant antibody, PBS Only (Capture) 83151-1-PBS (230329C5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TIMP1 Recombinant antibody, PBS Only (Detector) 83151-2-PBS (230329F11). 100 μ g. Concentration 1 mg/mL.

Alternative TIMP1 matched antibody pairs: MP00245-1, MP00245-3, MP50084-1, MP50084-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00245-2, TIMP1 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83151-1-PBS. Detection antibody: 83151-2-PBS. Standard: Eg0638. Range: 0.625-20 ng/mL