

For Research Use Only

TNFRSF11A Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00249-1

Capture Antibody Information

Catalog Number:
83253-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240106E7
Reactivity:
Human
GenBank:
NM_003839
Immunogen Catalog Number:
Ag29425

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 11a, NFKB activator
Gene ID:
8792

Detection Antibody Information

Catalog Number:
83253-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240106B10
Reactivity:
Human
GenBank:
NM_003839
Immunogen Catalog Number:
Ag29425

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 11a, NFKB activator
Gene ID:
8792

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00249-1 targets TNFRSF11A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFRSF11A Recombinant antibody, PBS Only (Capture) 83253-2-PBS (240106E7). 100 µg. Concentration 1 mg/ml.

Detection antibody: TNFRSF11A Recombinant antibody, PBS Only (Detector) 83253-1-PBS (240106B10). 100 µg. Concentration 1 mg/ml.

Alternative TNFRSF11A matched antibody pairs: MP00249-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00249-1, TNFRSF 11A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83253-2-PBS. Detection antibody: 83253-1-PBS. Standard: Ag29425. Range: 2.5-40 ng/mL.