

TPR Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00258-2**

Capture Antibody Information

Catalog Number:
83177-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230413C9
Reactivity:
Human
GenBank:
NM_003292

Conjugate:
Unconjugated
Full name:
translocated promoter region (to activated MET oncogene)
Gene ID:
7175

Detection Antibody Information

Catalog Number:
83177-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230413C3
Reactivity:
Human
GenBank:
NM_003292

Conjugate:
Unconjugated
Full name:
translocated promoter region (to activated MET oncogene)
Gene ID:
7175

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00258-2 targets TPR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TPR Recombinant antibody, PBS Only (Capture) 83177-4-PBS (230413C9). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TPR Recombinant antibody, PBS Only (Detector) 83177-3-PBS (230413C3). 100 μ g. Concentration 1 mg/mL.

Alternative TPR matched antibody pairs: MP00258-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00258-2, TPR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83177-4-PBS. Detection antibody: 83177-3-PBS. Standard: SY01141. Range: 0.312-5 ng/mL