

# DUS2L Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00261-1

## Capture Antibody Information

**Catalog Number:**  
83265-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240192B12  
**Reactivity:**  
Human  
**GenBank:**  
BC006527  
**Immunogen Catalog Number:**  
Ag23831

**Conjugate:**  
Unconjugated  
**Full name:**  
dihydrouridine synthase 2-like, SMM1 homolog (*S. cerevisiae*)  
**Gene ID:**  
54920

## Detection Antibody Information

**Catalog Number:**  
83265-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240192B1  
**Reactivity:**  
human  
**GenBank:**  
BC006527  
**Immunogen Catalog Number:**  
Ag23831

**Conjugate:**  
Unconjugated  
**Full name:**  
dihydrouridine synthase 2-like, SMM1 homolog (*S. cerevisiae*)  
**Gene ID:**  
54920

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.25-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00261-1 targets DUS2L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DUS2L Recombinant antibody, PBS Only (Capture) 83265-3-PBS (240192B12). 100 µg. Concentration 1 mg/ml.

Detection antibody: DUS2L Recombinant antibody, PBS Only (Detector) 83265-4-PBS (240192B1). 100 µg. Concentration 1 mg/ml.

Alternative DUS2L matched antibody pairs: MP00261-2, MP00261-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00261-1, DUS2L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83265-3-PBS. Detection antibody: 83265-4-PBS. Standard: Ag23831. Range: 1.25-40 ng/mL