

For Research Use Only



Rat IL-6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00264-3

Capture Antibody Information

Catalog Number:
83179-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240003G2

Reactivity:
Human

GenBank:
NM_012589.2

Conjugate:
Unconjugated

Full name:
interleukin 6

Gene ID:
24498

Detection Antibody Information

Catalog Number:
83179-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240003B1

Reactivity:
Human

GenBank:
NM_012589.2

Conjugate:
Unconjugated

Full name:
interleukin 6

Gene ID:
24498

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00264-3 targets IL-6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat IL-6 Recombinant antibody, PBS Only (Capture) 83179-2-PBS (240003G2). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Rat IL-6 Recombinant antibody, PBS Only (Capture/Detector) 83179-1-PBS (240003B1). 100 μ g. Concentration 1 mg/ml.

Alternative IL-6 matched antibody pairs: MP00264-1, MP00264-2, MP50032-2, MP50032-3, MP80002-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00264-3, Rat IL-6 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83179-2-PBS. Detection antibody: 83179-1-PBS. Standard: Eg0800. Range: 0.625-20 ng/mL