

BBS13 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00270-1

Capture Antibody Information

Catalog Number: 83181-1-PBS	Clone ID: 240036A2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: Meckel syndrome, type 1
Isotype: IgG	GenBank: BC010061	Gene ID: 54903
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9177	

Detection Antibody Information

Catalog Number: 83181-2-PBS	Clone ID: 240036E12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: Meckel syndrome, type 1
Isotype: IgG	GenBank: BC010061	Gene ID: 54903
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9177	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00270-1 targets BBS13 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BBS13 Recombinant antibody, PBS Only (Capture) 83181-1-PBS (240036A2). 100 µg. Concentration 1 mg/ml.

Detection antibody: BBS13 Recombinant antibody, PBS Only (Capture/Detector) 83181-2-PBS (240036E12). 100 µg. Concentration 1 mg/ml.

Alternative BBS13 matched antibody pairs: MP00270-2, MP00270-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Standard curve of MP00270-1, Human BBS13 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83181-1-PBS. Detection antibody: 83181-2-PBS. Standard: Ag9177. Range: 0.625-20 ng/mL