

BBS13 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00270-2

Capture Antibody Information

Catalog Number:
83181-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240036E12
Reactivity:
Human
GenBank:
BC010061
Immunogen Catalog Number:
Ag9177

Conjugate:
Unconjugated
Full name:
Meckel syndrome, type 1
Gene ID:
54903

Detection Antibody Information

Catalog Number:
83181-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240036E2
Reactivity:
Human
GenBank:
BC010061
Immunogen Catalog Number:
Ag9177

Conjugate:
Unconjugated
Full name:
Meckel syndrome, type 1
Gene ID:
54903

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00270-2 targets BBS13 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BBS13 Recombinant antibody, PBS Only (Capture/Detector) 83181-2-PBS (240036E12). 100 µg. Concentration 1 mg/ml.

Detection antibody: BBS13 Recombinant antibody, PBS Only (Detector) 83181-3-PBS (240036E2). 100 µg. Concentration 1 mg/ml.

Alternative BBS13 matched antibody pairs: MP00270-1, MP00270-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Standard curve of MP00270-2, Human BBS13 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83181-2-PBS. Detection antibody: 83181-3-PBS. Standard: Ag9177. Range: 0.625-20 ng/mL.