

GPR35 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00286-1

Capture Antibody Information

Catalog Number:
83197-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240032A11

Reactivity:
Human

GenBank:
NM_005301

Conjugate:
Unconjugated

Full name:
G protein-coupled receptor 35

Gene ID:
2859

Detection Antibody Information

Catalog Number:
83197-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240032B5

Reactivity:
Human

GenBank:
NM_005301

Conjugate:
Unconjugated

Full name:
G protein-coupled receptor 35

Gene ID:
2859

Applications

Tested Applications:
Cytometric bead array

Range:
78.125-5000 pg/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00286-1 targets GPR35 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GPR35 Recombinant antibody, PBS Only (Capture) 83197-1-PBS (240032A11). 100 µg. Concentration 1 mg/ml.

Detection antibody: GPR35 Recombinant antibody, PBS Only (Detector) 83197-3-PBS (240032B5). 100 µg. Concentration 1 mg/ml.

Alternative GPR35 matched antibody pairs: MP00286-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00286-1, GPR35 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83197-1-PBS. Detection antibody: 83197-3-PBS. Standard: SY01275. Range: 78.125-5000 pg/mL