

For Research Use Only



TFAP4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00288-3

Capture Antibody Information

Catalog Number:
83196-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240031C8
Reactivity:
Human
GenBank:
BC010576
Immunogen Catalog Number:
Ag2633

Conjugate:
Unconjugated
Full name:
transcription factor AP-4 (activating enhancer binding protein 4)
Gene ID:
7023

Detection Antibody Information

Catalog Number:
83196-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240031G6
Reactivity:
Human
GenBank:
BC010576
Immunogen Catalog Number:
Ag2633

Conjugate:
Unconjugated
Full name:
transcription factor AP-4 (activating enhancer binding protein 4)
Gene ID:
7023

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00288-3 targets TFAP4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TFAP4 Recombinant antibody, PBS Only (Capture/Detector) 83196-1-PBS (240031C8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TFAP4 Recombinant antibody, PBS Only (Detector) 83196-2-PBS (240031G6). 100 μ g. Concentration 1 mg/mL.

Alternative TFAP4 matched antibody pairs: MP00288-1, MP00288-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00288-3, TFAP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83196-1-PBS. Detection antibody: 83196-2-PBS. Standard: Ag2633. Range: 0.625-20 ng/mL