

For Research Use Only

GCNT1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00296-1

Capture Antibody Information

Catalog Number:
83121-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230447G12
Reactivity:
Human
GenBank:
BC074885
Immunogen Catalog Number:
Ag12264

Conjugate:
Unconjugated
Full name:
glucosaminyl (N-acetyl) transferase 1, core 2 (beta-1,6-N-acetylglucosaminyltransferase)
Gene ID:
2650

Detection Antibody Information

Catalog Number:
83121-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230447H7
Reactivity:
Human
GenBank:
BC074885
Immunogen Catalog Number:
Ag12264

Conjugate:
Unconjugated
Full name:
glucosaminyl (N-acetyl) transferase 1, core 2 (beta-1,6-N-acetylglucosaminyltransferase)
Gene ID:
2650

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00296-1 targets GCNT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GCNT1 Recombinant antibody, PBS Only (Capture) 83121-2-PBS (230447G12). 100 µg. Concentration 1 mg/ml.

Detection antibody: GCNT1 Recombinant antibody, PBS Only (Detector) 83121-3-PBS (230447H7). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00296-1, GCNT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83121-2-PBS. Detection antibody: 83121-3-PBS. Standard: Ag12264. Range: 0.625-20 ng/mL