

TGS1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00299-1

Capture Antibody Information

Catalog Number:
83140-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230427F6
Reactivity:
Human
GenBank:
BC011999
Immunogen Catalog Number:
Ag3683

Conjugate:
Unconjugated
Full name:
trimethylguanosine synthase homolog (*S. cerevisiae*)
Gene ID:
96764

Detection Antibody Information

Catalog Number:
83140-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230427G7
Reactivity:
Human
GenBank:
BC011999
Immunogen Catalog Number:
Ag3683

Conjugate:
Unconjugated
Full name:
trimethylguanosine synthase homolog (*S. cerevisiae*)
Gene ID:
96764

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00299-1 targets TGS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TGS1 Recombinant antibody, PBS Only (Capture) 83140-2-PBS (230427F6). 100 µg. Concentration 1 mg/ml.

Detection antibody: TGS1 Recombinant antibody, PBS Only (Detector) 83140-4-PBS (230427G7). 100 µg. Concentration 1 mg/ml.

Alternative TGS1 matched antibody pairs: MP00299-2, MP00299-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00299-1, TGS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83140-2-PBS. Detection antibody: 83140-4-PBS. Standard: Ag3683. Range: 0.625-20 ng/mL