

For Research Use Only

HMGB1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00311-1

Capture Antibody Information

Catalog Number: 82973-1-PBS	Clone ID: 230182G6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human, mouse, rat	Full name: high-mobility group box 1
Isotype: IgG	GenBank: BC003378	Gene ID: 3146
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1264	

Detection Antibody Information

Catalog Number: 82973-3-PBS	Clone ID: 230182C2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: high-mobility group box 1
Isotype: IgG	GenBank: BC003378	Gene ID: 3146
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1264	

Applications

Tested Applications: Sandwich ELISA	Range: 0.156-10 ng/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	---

Product Information

MP00311-1 targets HMGB1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.
Capture antibody: HMGB1 Recombinant antibody, PBS Only (Capture) 82973-1-PBS (230182G6). 100 µg. Concentration 1 mg/ml.

Detection antibody: HMGB1 Recombinant antibody, PBS Only (Detector) 82973-3-PBS (230182C2). 100 µg. Concentration 1 mg/ml.

Alternative HMGB1 matched antibody pairs: MP00311-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00311-1, HMGB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82973-1-PBS. Detection antibody: 82973-3-PBS. Standard: Ag1264. Range: 0.156-10 ng/mL