

For Research Use Only

GINS1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00317-1

Capture Antibody Information

Catalog Number:
83293-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240234B2
Reactivity:
Human
GenBank:
BC012542
Immunogen Catalog Number:
Ag8895

Conjugate:
Unconjugated
Full name:
GINS1 complex subunit 1 (Psf1 homolog)
Gene ID:
9837

Detection Antibody Information

Catalog Number:
83293-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240234F10
Reactivity:
Human
GenBank:
BC012542
Immunogen Catalog Number:
Ag8895

Conjugate:
Unconjugated
Full name:
GINS1 complex subunit 1 (Psf1 homolog)
Gene ID:
9837

Applications

Tested Applications:
Cytometric bead array

Range:
3.13-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00317-1 targets GINS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GINS1 Recombinant antibody, PBS Only (Capture) 83293-3-PBS (240234B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: GINS1 Recombinant antibody, PBS Only (Detector) 83293-4-PBS (240234F10). 100 µg. Concentration 1 mg/mL.

Alternative GINS1 matched antibody pairs: MP00317-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00317-1, GINS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83293-3-PBS. Detection antibody: 83293-4-PBS. Standard: Ag8895. Range: 3.13-100 ng/mL