

For Research Use Only

DDX47 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00325-1

Capture Antibody Information

Catalog Number:
83287-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240193E2
Reactivity:
Human
GenBank:
BC068009
Immunogen Catalog Number:
Ag11953

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-Asp) box polypeptide 47
Gene ID:
51202

Detection Antibody Information

Catalog Number:
83287-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240193A6
Reactivity:
Human
GenBank:
BC068009
Immunogen Catalog Number:
Ag11953

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-Asp) box polypeptide 47
Gene ID:
51202

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00325-1 targets DDX47 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDX47 Recombinant antibody, PBS Only (Capture) 83287-2-PBS (240193E2). 100 µg. Concentration 1 mg/ml.

Detection antibody: DDX47 Recombinant antibody, PBS Only (Detector) 83287-1-PBS (240193A6). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00325-1, DDX47 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83287-2-PBS. Detection antibody: 83287-1-PBS. Standard: Ag11953. Range: 0.625-80 ng/mL