

OX40L/TNFSF4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00327-2

Capture Antibody Information

Catalog Number:
82794-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240199F7
Reactivity:
Human
GenBank:
NM_003326.5

Conjugate:
Unconjugated
Full name:
tumor necrosis factor (ligand) superfamily, member 4
Gene ID:
7292

Detection Antibody Information

Catalog Number:
82794-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240199G9
Reactivity:
Human
GenBank:
NM_003326.5

Conjugate:
Unconjugated
Full name:
tumor necrosis factor (ligand) superfamily, member 4
Gene ID:
7292

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00327-2 targets OX40L/TNFSF4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OX40L/TNFSF4 Recombinant antibody, PBS Only (Capture) 82794-4-PBS (240199F7). 100 µg. Concentration 1 mg/ml.

Detection antibody: OX40L/TNFSF4 Recombinant antibody, PBS Only (Detector) 82794-3-PBS (240199G9). 100 µg. Concentration 1 mg/ml.

Alternative OX40L/TNFSF4 matched antibody pairs: MP00327-1, MP50147-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00327-2, OX40L/TNFSF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82794-4-PBS. Detection antibody: 82794-3-PBS. Standard: Eg0582. Range: 0.313-40ng/mL.