

For Research Use Only

# BCDIN3D Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00333-1

## Capture Antibody Information

<b>Catalog Number:</b> 83309-1-PBS	<b>Clone ID:</b> 240152A9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> BCDIN3 domain containing
<b>Isotype:</b> IgG	<b>GenBank:</b> BC053560	<b>Gene ID:</b> 144233
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag15613	

## Detection Antibody Information

<b>Catalog Number:</b> 83309-3-PBS	<b>Clone ID:</b> 240152C12	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> BCDIN3 domain containing
<b>Isotype:</b> IgG	<b>GenBank:</b> BC053560	<b>Gene ID:</b> 144233
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag15613	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.313-40 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP00333-1 targets BCDIN3D in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCDIN3D Recombinant antibody, PBS Only (Capture) 83309-1-PBS (240152A9). 100 µg. Concentration 1 mg/mL.

Detection antibody: BCDIN3D Recombinant antibody, PBS Only (Detector) 83309-3-PBS (240152C12). 100 µg. Concentration 1 mg/mL.

Alternative BCDIN3D matched antibody pairs: MP00333-2, MP00333-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00333-1, BCDIN3D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83309-1-PBS. Detection antibody: 83309-3-PBS. Standard: Ag15613. Range: 0.313-40 ng/mL.