

For Research Use Only

SMARCA4 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00334-2

Capture Antibody Information

Catalog Number:
83310-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240172D5
Reactivity:
Human
GenBank:
BC150298
Immunogen Catalog Number:
Ag16256

Conjugate:
Unconjugated
Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4
Gene ID:
6597

Detection Antibody Information

Catalog Number:
83310-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240172E7
Reactivity:
Human
GenBank:
BC150298
Immunogen Catalog Number:
Ag16256

Conjugate:
Unconjugated
Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4
Gene ID:
6597

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00334-2 targets SMARCA4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SMARCA4 Recombinant antibody, PBS Only (Capture) 83310-1-PBS (240172D5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SMARCA4 Recombinant antibody, PBS Only (Detector) 83310-5-PBS (240172E7). 100 μ g. Concentration 1 mg/mL.

Alternative SMARCA4 matched antibody pairs: MP00334-1, MP00334-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00334-2, SMARCA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83310-1-PBS. Detection antibody: 83310-5-PBS. Standard: Ag16256. Range: 0.625-20 ng/mL.