

# BRAP Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00338-1

## Capture Antibody Information

<b>Catalog Number:</b> 83314-2-PBS	<b>Clone ID:</b> 240099H3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> BRCA1 associated protein
<b>Isotype:</b> IgG	<b>GenBank:</b> BC136698	<b>Gene ID:</b> 8315
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag26728	

## Detection Antibody Information

<b>Catalog Number:</b> 83314-5-PBS	<b>Clone ID:</b> 240099G12	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> BRCA1 associated protein
<b>Isotype:</b> IgG	<b>GenBank:</b> BC136698	<b>Gene ID:</b> 8315
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag26728	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.313-40 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	---------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

## Product Information

MP00338-1 targets BRAP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BRAP Recombinant antibody, PBS Only (Capture) 83314-2-PBS (240099H3). 100 µg. Concentration 1 mg/mL.

Detection antibody: BRAP Recombinant antibody, PBS Only (Detector) 83314-5-PBS (240099G12). 100 µg. Concentration 1 mg/mL.

Alternative BRAP matched antibody pairs: MP00338-2, MP00338-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00338-1, BRAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83314-2-PBS. Detection antibody: 83314-5-PBS. Standard: Ag26728. Range: 0.313-40 ng/mL