

For Research Use Only

BRAP Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00338-2

Capture Antibody Information

Catalog Number: 83314-3-PBS	Clone ID: 240099A7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: BRCA1 associated protein
Isotype: IgG	GenBank: BC136698	Gene ID: 8315
Purification Method: Protein A purification	Immunogen Catalog Number: Ag26728	

Detection Antibody Information

Catalog Number: 83314-4-PBS	Clone ID: 240099A5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: BRCA1 associated protein
Isotype: IgG	GenBank: BC136698	Gene ID: 8315
Purification Method: Protein A purification	Immunogen Catalog Number: Ag26728	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00338-2 targets BRAP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BRAP Recombinant antibody, PBS Only (Capture) 83314-3-PBS (240099A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: BRAP Recombinant antibody, PBS Only (Detector) 83314-4-PBS (240099A5). 100 µg. Concentration 1 mg/mL.

Alternative BRAP matched antibody pairs: MP00338-1, MP00338-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00338-2, BRAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83314-3-PBS. Detection antibody: 83314-4-PBS. Standard: Ag26728. Range: 0.313-40 ng/mL