

For Research Use Only

SNAI1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00343-1

Capture Antibody Information

Catalog Number: 81584-2-PBS	Clone ID: 240164A5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: snail homolog 1 (Drosophila)
Isotype: IgG	GenBank: BC012910	Gene ID: 6615
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24248	

Detection Antibody Information

Catalog Number: 81584-3-PBS	Clone ID: 240164D2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: snail homolog 1 (Drosophila)
Isotype: IgG	GenBank: BC012910	Gene ID: 6615
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24248	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00343-1 targets SNAI1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SNAI1 Recombinant antibody, PBS Only (Capture) 81584-2-PBS (240164A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: SNAI1 Recombinant antibody, PBS Only (Detector) 81584-3-PBS (240164D2). 100 µg. Concentration 1 mg/ml.

Alternative SNAI1 matched antibody pairs: MP00343-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00343-1, SNAI1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81584-2-PBS. Detection antibody: 81584-3-PBS. Standard: Ag24248. Range: 0.313-20 ng/mL