

For Research Use Only

SAT1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00344-1

Capture Antibody Information

Catalog Number:
83319-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240157C6
Reactivity:
Human
GenBank:
BC008424
Immunogen Catalog Number:
Ag1154

Conjugate:
Unconjugated
Full name:
spermidine/spermine N1-acetyltransferase 1
Gene ID:
6303

Detection Antibody Information

Catalog Number:
83319-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240157H3
Reactivity:
Human
GenBank:
BC008424
Immunogen Catalog Number:
Ag1154

Conjugate:
Unconjugated
Full name:
spermidine/spermine N1-acetyltransferase 1
Gene ID:
6303

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00344-1 targets SAT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SAT1 Recombinant antibody, PBS Only (Capture) 83319-3-PBS (240157C6). 100 µg. Concentration 1 mg/ml.

Detection antibody: SAT1 Recombinant antibody, PBS Only (Capture/Detector) 83319-2-PBS (240157H3). 100 µg. Concentration 1 mg/ml.

Alternative SAT1 matched antibody pairs: MP00344-2, MP00344-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00344-1, SAT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83319-3-PBS. Detection antibody: 83319-2-PBS. Standard: Ag1154. Range: 0.625-80 ng/mL