

For Research Use Only

APIP Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00346-2

Capture Antibody Information

Catalog Number: 83320-4-PBS	Clone ID: 240189H7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: APAF1 interacting protein
Isotype: IgG	GenBank: BC009077	Gene ID: 51074
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2729	

Detection Antibody Information

Catalog Number: 83320-1-PBS	Clone ID: 240189A9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: APAF1 interacting protein
Isotype: IgG	GenBank: BC009077	Gene ID: 51074
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2729	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00346-2 targets APIP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: APIP Recombinant antibody, PBS Only (Capture) 83320-4-PBS (240189H7). 100 µg. Concentration 1 mg/ml.

Detection antibody: APIP Recombinant antibody, PBS Only (Detector) 83320-1-PBS (240189A9). 100 µg. Concentration 1 mg/ml.

Alternative APIP matched antibody pairs: MP00346-1, MP00346-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00346-2, AIP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83320-4-PBS. Detection antibody: 83320-1-PBS. Standard: Ag2729. Range: 0.313-40 ng/mL