

For Research Use Only

DNASE1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00356-1

Capture Antibody Information

Catalog Number: 83337-1-PBS	Clone ID: 240222A6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: deoxyribonuclease I
Isotype: IgG	GenBank: BC029437	Gene ID: 1773
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10117	

Detection Antibody Information

Catalog Number: 83337-5-PBS	Clone ID: 240222B11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: deoxyribonuclease I
Isotype: IgG	GenBank: BC029437	Gene ID: 1773
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10117	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00356-1 targets DNASE1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DNASE1 Recombinant antibody, PBS Only (Capture) 83337-1-PBS (240222A6). 100 µg. Concentration 1 mg/ml.

Detection antibody: DNASE1 Recombinant antibody, PBS Only (Detector) 83337-5-PBS (240222B11). 100 µg. Concentration 1 mg/ml.

Alternative DNASE1 matched antibody pairs: MP00356-2, MP00356-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00356-1, DNASE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83337-1-PBS. Detection antibody: 83337-5-PBS. Standard: Ag10117. Range: 0.313-40 ng/mL.