

For Research Use Only



# AGPAT2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00360-1

## Capture Antibody Information

Catalog Number:  
83349-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240153F11

Reactivity:  
Human

GenBank:  
BC019292

Immunogen Catalog Number:  
Ag34761

Conjugate:  
Unconjugated

Full name:  
1-acylglycerol-3-phosphate O-acyltransferase 2 (lysophosphatidic acid acyltransferase, beta)

Gene ID:  
10555

## Detection Antibody Information

Catalog Number:  
83349-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240153H1

Reactivity:  
Human

GenBank:  
BC019292

Immunogen Catalog Number:  
Ag34761

Conjugate:  
Unconjugated

Full name:  
1-acylglycerol-3-phosphate O-acyltransferase 2 (lysophosphatidic acid acyltransferase, beta)

Gene ID:  
10555

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00360-1 targets AGPAT2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AGPAT2 Recombinant antibody, PBS Only (Capture) 83349-2-PBS (240153F11). 100 µg. Concentration 1 mg/ml.

Detection antibody: AGPAT2 Recombinant antibody, PBS Only (Capture/Detector) 83349-4-PBS (240153H1). 100 µg. Concentration 1 mg/ml.

Alternative AGPAT2 matched antibody pairs: MP00360-2, MP00360-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00360-1, AGPAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83349-2-PBS. Detection antibody: 83349-4-PBS. Standard: Ag34761. Range: 0.625-80 ng/mL.