

CACHD1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00362-1

Capture Antibody Information

Catalog Number: 83350-2-PBS	Clone ID: 240197G8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: cache domain containing 1
Isotype: IgG	GenBank: BC039301	Gene ID: 57685
Purification Method: Protein A purification	Immunogen Catalog Number: Ag22138	

Detection Antibody Information

Catalog Number: 83350-1-PBS	Clone ID: 240197B2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: cache domain containing 1
Isotype: IgG	GenBank: BC039301	Gene ID: 57685
Purification Method: Protein A purification	Immunogen Catalog Number: Ag22138	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00362-1 targets CACHD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CACHD1 Recombinant antibody, PBS Only (Capture) 83350-2-PBS (240197G8). 100 µg. Concentration 1 mg/ml.

Detection antibody: CACHD1 Recombinant antibody, PBS Only (Detector) 83350-1-PBS (240197B2). 100 µg. Concentration 1 mg/ml.

Alternative CACHD1 matched antibody pairs: MP00362-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00362-1, CACHD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83350-2-PBS. Detection antibody: 83350-1-PBS. Standard: Ag22138. Range: 0.625-80 ng/mL.