

TCIRG1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00365-3

Capture Antibody Information

Catalog Number:
83351-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240166A4
Reactivity:
Human
GenBank:
BC018133
Immunogen Catalog Number:
Ag3378

Conjugate:
Unconjugated
Full name:
T-cell, immune regulator 1, ATPase, H⁺ transporting, lysosomal VO subunit A3
Gene ID:
10312

Detection Antibody Information

Catalog Number:
83351-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240166G12
Reactivity:
Human
GenBank:
BC018133
Immunogen Catalog Number:
Ag3378

Conjugate:
Unconjugated
Full name:
T-cell, immune regulator 1, ATPase, H⁺ transporting, lysosomal VO subunit A3
Gene ID:
10312

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00365-3 targets TCIRG1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TCIRG1 Recombinant antibody, PBS Only (Capture) 83351-1-PBS (240166A4). 100 µg. Concentration 1 mg/ml.

Detection antibody: TCIRG1 Recombinant antibody, PBS Only (Detector) 83351-2-PBS (240166G12). 100 µg. Concentration 1 mg/ml.

Alternative TCIRG1 matched antibody pairs: MP00365-1, MP00365-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00365-3, TCIRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83351-1-PBS. Detection antibody: 83351-2-PBS. Standard: Ag3378. Range: 0.313-40 ng/mL