

# KIAA1199 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00373-1

## Capture Antibody Information

**Catalog Number:**  
83373-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240372C3  
**Reactivity:**  
Human  
**GenBank:**  
BC020256  
**Immunogen Catalog Number:**  
Ag15527

**Conjugate:**  
Unconjugated  
**Full name:**  
KIAA1199  
**Gene ID:**  
57214

## Detection Antibody Information

**Catalog Number:**  
83373-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240372B7  
**Reactivity:**  
Human  
**GenBank:**  
BC020256  
**Immunogen Catalog Number:**  
Ag15527

**Conjugate:**  
Unconjugated  
**Full name:**  
KIAA1199  
**Gene ID:**  
57214

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00373-1 targets KIAA1199 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KIAA1199 Recombinant antibody, PBS Only (Capture) 83373-3-PBS (240372C3). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: KIAA1199 Recombinant antibody, PBS Only (Capture/Detector) 83373-1-PBS (240372B7). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative KIAA1199 matched antibody pairs: MP00373-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00373-1, KIAA1199 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83373-3-PBS. Detection antibody: 83373-1-PBS. Standard: Ag15527. Range: 0.625-80 ng/mL.