

For Research Use Only

RAT Ccl2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00390-1

Capture Antibody Information

Catalog Number:
83370-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240341E10

Reactivity:
Rat

GenBank:
NM_031530

Conjugate:
Unconjugated

Full name:
chemokine (C-C motif) ligand 2

Gene ID:
24770

Detection Antibody Information

Catalog Number:
83370-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240341D11

Reactivity:
Rat

GenBank:
NM_031530

Conjugate:
Unconjugated

Full name:
chemokine (C-C motif) ligand 2

Gene ID:
24770

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00390-1 targets Ccl2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat MCP-1 Recombinant antibody, PBS Only (Capture) 83370-3-PBS (240341E10). 100 µg. Concentration 1 mg/ml.

Detection antibody: Rat MCP-1 Recombinant antibody, PBS Only (Detector) 83370-1-PBS (240341D11). 100 µg. Concentration 1 mg/ml.

Alternative Ccl2 matched antibody pairs: MP00168-1, MP00168-2, MP00168-3, MP00390-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00390-1, RAT MCP-1/CCL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83370-3-PBS. Detection antibody: 83370-1-PBS. Standard: Eg0397. Range: 0.313-20 ng/mL.