

For Research Use Only

DUSP6 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00391-4

Capture Antibody Information

Catalog Number:
83356-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240295E12

Reactivity:
Human

GenBank:
BC003562

Immunogen Catalog Number:
Ag23474

Conjugate:
Unconjugated

Full name:
dual specificity phosphatase 6

Gene ID:
1848

Detection Antibody Information

Catalog Number:
83356-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240295G11

Reactivity:
Human

GenBank:
BC003562

Immunogen Catalog Number:
Ag23474

Conjugate:
Unconjugated

Full name:
dual specificity phosphatase 6

Gene ID:
1848

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00391-4 targets DUSP6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DUSP6 Recombinant antibody, PBS Only (Capture) 83356-3-PBS (240295E12). 100 μ g. Concentration 1 mg/mL.

Detection antibody: DUSP6 Recombinant antibody, PBS Only (Detector) 83356-4-PBS (240295G11). 100 μ g. Concentration 1 mg/mL.

Alternative DUSP6 matched antibody pairs: MP00391-1, MP00391-2, MP00391-3, MP00391-5, MP00391-6

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00391-4, DUSP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83356-3-PBS. Detection antibody: 83356-4-PBS. Standard: Ag23474. Range: 0.625-80 ng/mL