

TLN1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00405-3

Capture Antibody Information

Catalog Number:
82856-16-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240294C11
Reactivity:
Human
GenBank:
BC042923
Immunogen Catalog Number:
Ag5367

Conjugate:
Unconjugated
Full name:
talin 1
Gene ID:
7094

Detection Antibody Information

Catalog Number:
82856-15-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240294F4
Reactivity:
Human
GenBank:
BC042923
Immunogen Catalog Number:
Ag5367

Conjugate:
Unconjugated
Full name:
talin 1
Gene ID:
7094

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00405-3 targets TLN1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TLN1 Recombinant antibody, PBS Only (Capture) 82856-16-PBS (240294C11). 100 µg. Concentration 1 mg/ml.

Detection antibody: TLN1 Recombinant antibody, PBS Only (Detector) 82856-15-PBS (240294F4). 100 µg. Concentration 1 mg/ml.

Alternative TLN1 matched antibody pairs: MP00405-1, MP00405-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00405-3, Talin-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82856-16-PBS. Detection antibody: 82856-15-PBS. Standard: Ag5367. Range: 1.56-100 ng/mL