

FGFR3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00407-2

Capture Antibody Information

Catalog Number:
83701-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240758G6

Reactivity:
Human

GenBank:
NM_001163213.1

Conjugate:
Unconjugated

Full name:
fibroblast growth factor receptor 3

Gene ID:
2261

Detection Antibody Information

Catalog Number:
83701-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240758B1

Reactivity:
Human

GenBank:
NM_001163213.1

Conjugate:
Unconjugated

Full name:
fibroblast growth factor receptor 3

Gene ID:
2261

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00407-2 targets FGFR3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FGFR3 Recombinant antibody, PBS Only (Capture) 83701-2-PBS (240758G6). 100 µg. Concentration 1 mg/ml.

Detection antibody: FGFR3 Recombinant antibody, PBS Only (Detector) 83701-1-PBS (240758B1). 100 µg. Concentration 1 mg/ml.

Alternative FGFR3 matched antibody pairs: MP00407-1, MP00407-3, MP00407-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00407-2, FGFR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83701-2-PBS. Detection antibody: 83701-1-PBS. Standard: Eg1041. Range: 0.313-40 ng/mL