

For Research Use Only

Mouse Il22 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00410-2

Capture Antibody Information

Catalog Number:
83405-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240251A12

Reactivity:
Mouse

GenBank:
NM_016971

Conjugate:
Unconjugated

Full name:
interleukin 22

Gene ID:
50929

Detection Antibody Information

Catalog Number:
83405-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240251F9

Reactivity:
Mouse

GenBank:
NM_016971

Conjugate:
Unconjugated

Full name:
interleukin 22

Gene ID:
50929

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00410-2 targets IL22 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse IL22 Recombinant antibody, PBS Only (Capture) 83405-1-PBS (240251A12). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Mouse IL22 Recombinant antibody, PBS Only (Capture/Detector) 83405-2-PBS (240251F9). 100 μ g. Concentration 1 mg/ml.

Alternative IL22 matched antibody pairs: MP00410-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00410-2, MOUSE IL-22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83405-1-PBS. Detection antibody: 83405-2-PBS. Standard: Eg0405. Range: 0.625-40 ng/mL.