

FEZ2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00422-1

Capture Antibody Information

Catalog Number:
83407-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240248C9
Reactivity:
Human
GenBank:
BC110386
Immunogen Catalog Number:
Ag15370

Conjugate:
Unconjugated
Full name:
fasciculation and elongation protein zeta 2 (zygin II)
Gene ID:
9637

Detection Antibody Information

Catalog Number:
83407-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240248G8
Reactivity:
Human
GenBank:
BC110386
Immunogen Catalog Number:
Ag15370

Conjugate:
Unconjugated
Full name:
fasciculation and elongation protein zeta 2 (zygin II)
Gene ID:
9637

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00422-1 targets FEZ2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FEZ2 Recombinant antibody, PBS Only (Capture) 83407-2-PBS (240248C9). 100 µg. Concentration 1 mg/ml.

Detection antibody: FEZ2 Recombinant antibody, PBS Only (Detector) 83407-5-PBS (240248G8). 100 µg. Concentration 1 mg/ml.

Alternative FEZ2 matched antibody pairs: MP00422-2, MP00422-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00422-1, FEZ2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83407-2-PBS. Detection antibody: 83407-5-PBS. Standard: Ag15370. Range: 0.313-40 ng/mL