

For Research Use Only

GNB3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00443-3

Capture Antibody Information

Catalog Number:
83452-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240281B5
Reactivity:
Human
GenBank:
BC000115
Immunogen Catalog Number:
Ag0127

Conjugate:
Unconjugated
Full name:
guanine nucleotide binding protein (G protein), beta polypeptide 3
Gene ID:
2784

Detection Antibody Information

Catalog Number:
83452-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240281G11
Reactivity:
Human
GenBank:
BC000115
Immunogen Catalog Number:
Ag0127

Conjugate:
Unconjugated
Full name:
guanine nucleotide binding protein (G protein), beta polypeptide 3
Gene ID:
2784

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00443-3 targets GNB3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GNB3 Recombinant antibody, PBS Only (Capture/Detector) 83452-1-PBS (240281B5). 100 µg. Concentration 1 mg/mL.

Detection antibody: GNB3 Recombinant antibody, PBS Only (Detector) 83452-2-PBS (240281G11). 100 µg. Concentration 1 mg/mL.

Alternative GNB3 matched antibody pairs: MP00443-1, MP00443-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00443-3, GNB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83452-1-PBS. Detection antibody: 83452-2-PBS. Standard: Ag0127. Range: 0.78-100 ng/mL