

PLXDC2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00445-3

Capture Antibody Information

Catalog Number:
83449-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240174F1
Reactivity:
Human
GenBank:
BC012885
Immunogen Catalog Number:
Ag2940

Conjugate:
Unconjugated
Full name:
plexin domain containing 2
Gene ID:
84898

Detection Antibody Information

Catalog Number:
83449-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240174A3
Reactivity:
Human
GenBank:
BC012885
Immunogen Catalog Number:
Ag2940

Conjugate:
Unconjugated
Full name:
plexin domain containing 2
Gene ID:
84898

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00445-3 targets PLXDC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PLXDC2 Recombinant antibody, PBS Only (Capture/Detector) 83449-2-PBS (240174F1). 100 µg. Concentration 1 mg/ml.

Detection antibody: PLXDC2 Recombinant antibody, PBS Only (Detector) 83449-1-PBS (240174A3). 100 µg. Concentration 1 mg/ml.

Alternative PLXDC2 matched antibody pairs: MP00445-1, MP00445-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00445-3, PLXDC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83449-2-PBS. Detection antibody: 83449-1-PBS. Standard: Ag2940. Range: 0.313-40 ng/mL