

For Research Use Only

# HGF Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00454-2

## Capture Antibody Information

**Catalog Number:**  
83453-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240383C7  
**Reactivity:**  
Human  
**GenBank:**  
BC130286  
**Immunogen Catalog Number:**  
Ag25478

**Conjugate:**  
Unconjugated  
**Full name:**  
hepatocyte growth factor  
(hepapoietin A; scatter factor)  
**Gene ID:**  
3082

## Detection Antibody Information

**Catalog Number:**  
83453-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240383E11  
**Reactivity:**  
Human  
**GenBank:**  
BC130286  
**Immunogen Catalog Number:**  
Ag25478

**Conjugate:**  
Unconjugated  
**Full name:**  
hepatocyte growth factor  
(hepapoietin A; scatter factor)  
**Gene ID:**  
3082

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00454-2 targets HGF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HGF Recombinant antibody, PBS Only (Capture) 83453-3-PBS (240383C7). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: HGF Recombinant antibody, PBS Only (Capture/Detector) 83453-2-PBS (240383E11). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative HGF matched antibody pairs: MP00454-1, MP50034-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00454-2, HGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83453-3-PBS. Detection antibody: 83453-2-PBS. Standard: Ag25478. Range: 0.313-40 ng/mL