

For Research Use Only

C2orf7 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00468-1

Capture Antibody Information

Catalog Number:
83478-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240494C7

Reactivity:
Human

GenBank:
BC005069

Immunogen Catalog Number:
Ag14678

Conjugate:
Unconjugated

Full name:
chromosome 2 open reading frame 7

Gene ID:
84279

Detection Antibody Information

Catalog Number:
83478-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240494G2

Reactivity:
Human

GenBank:
BC005069

Immunogen Catalog Number:
Ag14678

Conjugate:
Unconjugated

Full name:
chromosome 2 open reading frame 7

Gene ID:
84279

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00468-1 targets C2orf7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C2orf7 Recombinant antibody, PBS Only (Capture) 83478-3-PBS (240494C7). 100 µg. Concentration 1 mg/ml.

Detection antibody: C2orf7 Recombinant antibody, PBS Only (Detector) 83478-2-PBS (240494G2). 100 µg. Concentration 1 mg/ml.

Alternative C2orf7 matched antibody pairs: MP00468-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00468-1, C2orf7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83478-3-PBS. Detection antibody: 83478-2-PBS. Standard: Ag14678. Range: 1.25-80 ng/mL