

TFF1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00469-1

Capture Antibody Information

Catalog Number:
83467-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240424F5

Reactivity:
Human

GenBank:
BC032811

Immunogen Catalog Number:
Ag4541

Conjugate:
Unconjugated

Full name:
trefoil factor 1

Gene ID:
7031

Detection Antibody Information

Catalog Number:
83467-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240424E6

Reactivity:
Human

GenBank:
BC032811

Immunogen Catalog Number:
Ag4541

Conjugate:
Unconjugated

Full name:
trefoil factor 1

Gene ID:
7031

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00469-1 targets TFF1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TFF1 Recombinant antibody, PBS Only (Capture) 83467-2-PBS (240424F5). 100 µg. Concentration 1 mg/ml.

Detection antibody: TFF1 Recombinant antibody, PBS Only (Detector) 83467-1-PBS (240424E6). 100 µg. Concentration 1 mg/ml.

Alternative TFF1 matched antibody pairs: MP50027-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00469-1, TFF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83467-2-PBS. Detection antibody: 83467-1-PBS. Standard: Ag4541. Range: 0.313-40 ng/mL